



**Product: Enema Set**  
**REF: 61001 – Irrigation Set**

**Specifications:**

<b>Product:</b>	
Bag dimensions	17 x 35 cm
Bag material	DEHP free PVC Softener: DINP
Tube length	120 cm
Tube dimensions (outer/inner dia)	Ø9.0 mm x Ø 6.6 mm
Tube material	DEHP free PVC
Tube tip	DEHP free PVC
Tube tip – dimensions (outer/inner dia & length)	Ø10.0 x Ø 7.4 x 150 mm
Valve	POM

<b>Packaging:</b>	
Pouch – single pack	LDPE - 17 cm x 27 cm
Bag (for 100 single packed products)	LDPE - 80 cm x 100 cm
Transport carton	100 single packed products

<b>IFU:</b>	
Instructions on the retail label	Yes

<b>Regulatory:</b>	
MDR – Class I (non-sterile)	DoC

<b>Shelf life:</b>	
From date of production	5 years